

# Manual of k-t FOCUSS with ME/MC

Revised: Nov. 06, 2012

Contact: [jong.ye@kaist.ac.kr](mailto:jong.ye@kaist.ac.kr) (Prof. Jong Chul Ye)

Homepage: <http://bisp.kaist.ac.kr>

This program is with respect to the *k-t FOCUSS with ME/MC*.

We divided it into two categories of Cartesian and radial formats, respectively.

Each format includes three main functions of measurement generation with down-sampling rate, *k-t FOCUSS* and *k-t FOCUSS with ME/MC*.

Execution files are placed in folder '*bin*', and example data is placed in folder '*data*'.

Here are more details.

A. Cartesian format (you can just execute Exe1=> Exe3)

1. '*Exe1\_Measurement\_Gen.m*'

- 1) Open the file
- 2) Define down-sampling rate: *DsRate*
- 3) Select reference image type for ME/MC calculation
  - a. A fully sampled frame: '*FullSingleFrame*'
  - b. Temporal average of diastole frames: '*DiastoleFrames*'
- 4) Define the number of central phase encoding lines: '*num\_low\_phase*'
- 5) In case 'a', choose one frame as a reference image
- 6) In case 'b', define diastole frames manually to generate a reference image
- 7) Execute the file

2. '*Exe2\_ktFOCUSS.m*'

- 1) Open the file
- 2) Define functions: *A*, *AT*

- 3) Define parameters for  $k$ - $t$  FOCUSS
- 4) Execute the file
3. 'Exe3\_ktFOCUSSwithMEMC.m'
- 1) Open the file
- 2) Define functions:  $A$ ,  $AT$
- 3) Define parameters for  $k$ - $t$  FOCUSS
- 4) Define parameters for ME/MC:  $px$ ,  $py$ ,  $ws$
- 5) Execute the file

#### B. Radial format (you can just execute Exe1=> Exe3)

Our radial data has three dimensions:  $(r, \phi, t)$ .

We first calculate the re-binning from radial to the Cartesian domain, which is equivalent to the 2D-fourier transform domain.

Since the re-binning process is calculated in the shifted domain, low frequency of the original data should be centered (see functions: 'fftshift', 'ifftshift').

When you define functions of 'A' and 'AT', shift operators should be included properly (see our example code).

##### 1. 'Exe1\_Measurement\_Gen.m'

- 1) Open the file
- 2) Define down-sampling rate:  $DsRate$
- 3) Select reference image type for ME/MC calculation
  - a. A fully sampled frame: '**FullSingleFrame**'
  - b. Temporal average of diastole frames: '**DiastoleFrames**'
- 4) In case 'a', choose one frame as a reference image
- 5) In case 'b', define diastole frames manually to generate a reference image
- 6) Execute the file
2. 'Exe2\_ktFOCUSS.m'
- 1) Open the file

- 2) Define functions:  $A, AT$
  - 3) Define parameters for  $k$ - $t$  FOCUS
  - 4) Execute the file
3. 'Exe3\_ktFOCUSwithMEMC.m'
- 1) Open the file
  - 2) Define functions:  $A, AT$
  - 3) Define parameters for  $k$ - $t$  FOCUS
  - 4) Define parameters for ME/MC:  $px, py, ws$
  - 5) Execute the file